





Solutions for a nanoscale world.™



di EnviroScope


Atomic Force Microscope (AFM) System



- **Advanced Environmental Control**
- **Flexibility, Ease of Use, and Performance**
- **Full Range of Imaging Modes**

di EnviroScope

Proven Environmental Control for Co



The EnviroScope™ atomic force microscope (AFM) combines enhanced environmental controls within a sealed hermetic sample chamber, and the world's most advanced scanning technology to deliver greater application flexibility to research and industrial laboratories looking for a competitive edge. The system allows observation of sample reactions in a variety of complex environmental conditions, including high vacuum, liquid, and gas purge and exchange, as well as high temperature and fluid heating. Based on the industry-standard Dimension™ AFM head, the EnviroScope scanner provides proven scanning probe performance and reliability while performing an extensive range of imaging modes, from contact AFM and magnetic force microscopy to patented techniques like TappingMode™ AFM and PhasImaging™.


- ▶ Hermetically sealed sample chamber
 - **Permits scanning in high-vacuum, gas, liquid, and temperature-controlled environments**
 - **Enables comprehensive sample reaction studies**
- ▶ Modular stage design
 - **Provides easy sample setup**
 - **Enables application specific sample preparation**
- ▶ Dimension AFM scanning head
 - **Proven performance over a wide range of environmental conditions**
 - **Utilizes patented TrakScan™ technology for spot-on laser-beam reflection**
- ▶ Specialized electrochemistry module
 - **Temperature control**
 - **Atmospheric control above fluid cells**
 - **Prevents reactivity of the fluid and sample with the gaseous environment**

Flexibility over Research Environments

The EnviroScope continues Veeco's tradition of outpacing conventional systems by imaging a sample as it reacts with environmental conditions enabling more comprehensive observation. The EnviroScope is specifically designed to be a no-fuss, high-performance tool for research in material sciences, electrochemistry, polymer technologies, life sciences, and other applications.

The sophisticated, sample chamber is hermetically sealed. It can be purged with inert gas to slow or inhibit reactive processes and reduce moisture effects. The easy-to-use modular sample stages and an exclusive cantilever replacement tool enable users to vary their environmental experiments. For example, the optional fluid/heater stage allows fluid to surround the sample for observation of growth, etching, or corrosion at ambient and elevated temperatures.

For temperature-sensitive applications, the EnviroScope provides a high-temperature stage for observing reactions or crystallization in various temperature ranges over time. Software specially designed for the sample heating system allows programmable image acquisition linked to temperature for multiple-image thermal experiments.



EnviroScope modular stages. A) standard, B) fluid/heater, C) high temperature heater.

Turnkey Options for Customized Research

The EnviroScope also has an optional all-in-one turnkey vacuum system. This system utilizes a turbo pump that is mounted directly below the custom isolation table and is properly isolated from the AFM. It also has a roughing pump that can be remotely positioned to ensure no interference during imaging.

The electrochemistry (EC) option includes a fluid cell. The chamber allows gas flow above the cell to control sample evaporation and chemistry. Studies of interest include corrosion, etching, electroplating, electroadsorption, and electrodesorption. This option can accommodate samples up to 1 inch and provides true reference electrodes in the EC cell.

Complex AFM Studies



Turnkey vacuum system.

Unmatched Scanning Technology

The EnviroScope scanner is based on the industry-leading Dimension AFM head. A conformal coating has been added to this head to provide protection from environmental changes in the chamber. In addition, a built-in temperature sensor monitors the piezo temperature and provides a warning and shut-off if the piezo exceeds safe temperature levels. This scanner incorporates patented TrakScan technology to ensure that the laser beam reflects from the same spot on the AFM cantilever during raster scanning.

In addition, the EnviroScope offers a choice of high-performance NanoScope controllers. The NanoScope IIIa controller combines advanced analog and digital circuit designs with premium software and hardware to precisely control the SPM. The superior performance and utility of the NanoScope IIIa controller has led to more publications than all other SPM controllers combined.



The NanoScope IV controller incorporates all of the same features, as well as up to ten-times-faster scanning, increased functionality, bandwidth, flexibility, and expandability. The NanoScope IV incorporates Quadrex™ technology with lock-in detection and advanced signal routing to enable the patented Phasemaging technique. Particularly applicable for advanced studies of mechanical properties, Quadrex allows the detection of variations in adhesion, viscoelasticity, and other properties by mapping the phase of the cantilever oscillation during a scan. Using either controller, the EnviroScope can scan from 90 microns X-Y and 4.8 microns in Z with full, 16-bit resolution on all scan waveforms and on each axis.

Complete Range of Imaging Modes


The EnviroScope combines its environmental capacities with Veeco's unmatched list of imaging techniques for the detailed measurement of the mechanical properties and magnetic/electrical field of the surface.

- Torsional Resonance Mode (TRmode™)
- TappingMode
- Contact AFM
- MFM
- LFM
- EFM
- Phasemaging
- Surface Potential Microscopy
- Force Distance/Force Volume Measurements

Check with factory for compatibilities.




Poly-sbs block co-polymer. 26 Celsius, 760 Torr (left). 130 Celsius, 2×10^{-5} Torr (right). 1 μm scans, 5nm Z range.




Cover image: DNA in fluid. 2 μ m scan imaged with EnviroScope AFM.


MFM image of video tape in air at ambient (a), and in vacuum at 10⁻⁵ torr (b) showing higher amplitude variations and resolution. 5 μ m scans.



a.



b.



c.

SPP Polymer imaged with EnviroScope high-temperature heating stage at ambient (a), 60°C (b) and 140°C (c) temperatures. 5 μ m scans.

EnviroScope Specifications

Sample stage range: – 6mm X-Y, 14mm Z sample movement range while maintaining vacuum

Sample size: – 30mm X-Y; 12mm Z

Sample leveling: – Automatic mechanical

Temperature range while imaging:

- 18.5°C in ambient environment
- Ambient to 275°C in vacuum (to 300°C nominal at the sample)
- Ambient to 60°C in fluid (including EC)
- Temperature stability $\pm 1^\circ\text{C}$

Scan size: – 90 μ m X-Y; 5 μ m Z

Noise level:

- <0.05nm -RMS at ambient pressure
- Laboratory noise limit for best performance: 62 dBc

Linearity: – Software-corrected

Optics:

- Resolution on integrated top-view video microscope 3-4 μ m
- Fixed ~0.5mm field of view at the sample

Ports: – Gas purging, vacuum, or interface plate for user customization

Electronics controller: – Digital Instruments NanoScope IIIa with Quadrex or NanoScope IV

Vacuum level: – 10⁻⁵ Torr (requires optional purchase)

Note: Performance specifications are typical and subject to change without notice.

LASER RADIATION
Do not stare into the beam or view directly with optical instruments
class 3R laser product
1.6mW Max. @ 650-695nm.

2002 HCS, LLC
Reorder No. 492-008-003



112 Robin Hill Road, Santa Barbara, CA 93117
805-967-1400 • 1-888-24-VEECO
www.veeco.com

B49, Rev A1, 2/15/05

© 2005 Veeco Instruments Inc. All rights reserved.

EnviroScope, Dimension, PhaseImaging, TrakScan, TRmode, and Quadrex are trademarks of Veeco Instruments Inc.

NanoScope is a registered trademark of Veeco Instruments Inc.

Worldwide Customer Support from the Industry Leader

Veeco Instruments Inc. provides solutions for nanoscale applications in the worldwide semiconductor, data storage, telecommunications/wireless and scientific research markets. Our Metrology products are used to measure at the nanoscale and our Process Equipment tools help create nanoscale devices. Veeco's manufacturing and engineering facilities are located in New York, New Jersey, California, Colorado, Arizona and Minnesota. Global sales and service offices are located throughout the United States, Europe, Japan and Asia Pacific.